

## Directions for use:

**Pro-Zyme™** Low Foam Enzymatic Detergent is specially formulated for automatic washers where ultra low foaming detergents are desired.

Add 1/8 to 1 fl. oz. of **Pro-Zyme™** Low Foam Enzymatic Detergent per gallon of water. Optimal temperature range is 70° F-140° F.

As a manual presoak: Soak instruments in diluted **Pro-Zyme™** Low Foam Enzymatic Detergent solution for a minimum of 5 minutes.

As an ultrasonic cleaner or presoak: Add diluted **Pro-Zyme™** Low Foam Enzymatic Detergent solution to ultrasonic unit. Immerse instruments in enzyme cleaning solution. Operate for a minimum of 5 minutes.

As an automatic washer detergent and/or presoak: Program machine to deliver 1/8 to 1 fl. oz. of **Pro-Zyme™** Low Foam Enzymatic Detergent per gallon of water for both pre-wash and regular wash operations. Allow at least 5 minutes cycle time.