

Directions for use

Reconditioning of Stainless Steel Surgical Instruments: Dilute one ounce of **Restore Sheen™** Rust & Stain Remover to 1 gallon of water in a stainless steel tray. Open clean instruments and completely immerse in diluted solution for 20 minutes or for 15 minutes in a fully heated ultrasonic washer. Severely stained instruments may require up to 30 minutes of soak time. Remove and rinse with water. If needed, treat again by scrubbing with **Restore Sheen™** Rust & Stain Remover, then rinsing thoroughly with water.

Severely stained instruments may require a second treatment.

Dilute one ounce of **Restore Sheen™** Rust & Stain Remover to 1 gallon of water. Brush stains and rust away with stiff brush and hot solution. Rinse completely following wash procedures.

DANGER - HAZARDS: May cause fire or explosion; strong oxidizer. May be corrosive to metals. Causes severe skin burns and eye damage. Causes serious eye damage. Harmful to aquatic life with long lasting effects.



PRECAUTIONS: Keep away from heat/sparks/open flames/hot surfaces. - No smoking. Keep away from combustible materials. Take any precaution to avoid mixing with combustibles. Keep only in original container.

P264 - Wash hands thoroughly after handling. Avoid release to the environment. Wear protective gloves/protective clothing/eye protection/face protection. Wear fire/flame resistant/retardant clothing

RESPONSE: IF SWALLOWED: Rinse mouth. Do not induce vomiting. Call a POISON CENTER or doctor/physician. **IF ON SKIN:** Wash with plenty of soap and water. If skin irritation occurs, get medical advice/attention. Take off contaminated clothing and wash before re-use. **IF INHALED:** Remove victim to fresh air and keep at rest in a position comfortable for breathing. **IF IN EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Call a poison control center or physician for treatment advice. Have product container or label with you when calling poison control center or physician. If eye irritation persists, get medical advice/attention.

In case of major fire and large quantities: evacuate area and fight fire remotely due to the risk of explosion. Absorb spillage to prevent material damage.

STORAGE & DISPOSAL:

Store locked up. Unused product is not a RCRA Hazardous waste. However, contaminated product and wastes may be RCRA hazardous. Users are advised to determine the appropriate disposal method based on local, state and federal regulations and comply with those regulations.